

Anti-CD5 Antibody (1Q411)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Rat
Conjugation:	Unconjugated
Clone:	1Q411
Purification:	Protein A

Applications

Verified Activity:	<p>1. Flow cytometric analysis of Rat CD5 expression on SD rat splenocytes. Cells were stained with purified anti-Rat CD5, then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.</p> <p>2. Immunofluorescence staining of ratCD5 in rat spleen cells. Cells were fixed with 4% PFA, blocked with 10% serum, and incubated with mouse anti-ratCD5 monoclonal antibody (dilution ratio 1:60) at 4°C overnight. Then cells were stained with the Alexa Fluor®488-conjugated Goat Anti-mouse IgG secondary antibody (green). Positive staining was localized to Cell membrane.</p>
Application:	FCM,ICC/IF
Recommended	ICC-IF: 1:20-1:100; FCM: 1:25-1:100

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Rat CD5 / Cluster of Differentiation 5 Protein (TMPY-03169)
Antigen Species:	Rat
Synonyms:	LEU1;T1;CD5 molecule
Biology Area:	ITIM/ITAM Immunoreceptors and Related Molecules

Research Background

The cluster of differentiation (CD) system is commonly used as cell markers in Immunophenotyping. Different kinds of cells in the immune system can be identified through the surface CD molecules associating with the immune function of the cell. There are more than 320 CD unique clusters and subclusters have been identified. Some of the CD molecules serve as receptors or ligands important to the cell through initiating a signal cascade which then alter the behavior of the cell. Some CD proteins do not take part in cell signal process but have other functions such as cell adhesion. CD5 is a member of the CD system. CD5 was found to be widely distributed in T-cells and B1 cells which is a subset of IgM-secreting B cells. CD5 also was found expressed in small lymphocytic lymphoma, hairy cell leukaemia and mantle cell lymphoma cells. CD5 serves to weaken the activating stimulus from the BCR so that the B1 cells can only reflect to the very strong stimuli but not the normal tissue proteins.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481