

## Anti-p53 Antibody (3K688)

### Product Details

Ig Type:	Rabbit IgG
Reactivity:	Cynomolgus, Human
Conjugation:	Unconjugated
Clone:	3K688
Purification:	Protein A

### Applications

Verified Activity:	<p>1. cyno p53 was immunoprecipitated using:</p> <ul style="list-style-type: none"><li>-Lane A:0.5 mg A431 Whole Cell Lysate</li><li>-0.5 µL anti-cyno p53 rabbit monoclonal antibody and 60 µg of Immunomagnetic beads Protein G.</li></ul> <p>-Primary antibody:</p> <ul style="list-style-type: none"><li>-Anti-cyno p53 rabbit monoclonal antibody, at 1:500 dilution.</li></ul> <p>-Secondary antibody:</p> <ul style="list-style-type: none"><li>-Clean-Blot<sup>®</sup> IP Detection Reagent (HRP) at 1:500 dilution.</li><li>-Developed using the DAB staining technique.</li><li>-Performed under reducing conditions.</li></ul> <p>-Predicted band size: 55 kDa.</p> <ul style="list-style-type: none"><li>-Observed band size: 55 kDa.</li></ul> <p>2. Anti-p53 rabbit monoclonal antibody at 1:500 dilution.</p> <ul style="list-style-type: none"><li>-Lane A: A431 Whole Cell lysate.</li><li>-Lysates/proteins at 30 µg per lane.</li></ul> <p>-Secondary</p> <ul style="list-style-type: none"><li>-Goat Anti-Rabbit IgG H&amp;L (Dylight800) at 1/10000 dilution.</li><li>-Developed using the Odyssey technique.</li><li>-Performed under reducing conditions.</li><li>-Predicted band size:53 kDa.</li><li>-Observed band size:53 kDa</li></ul>
Application:	ELISA,IP,WB
Recommended	WB: 1:500-1:1000; ELISA: 1:25000-1:50000; IP: 0.2-1 µL/mg of lysate

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Cynomolgus p53 / TP53 protein (TMPY-02228)
Antigen Species:	Cynomolgus
Synonyms:	tumor protein p53;TP53;p53
Biology Area:	Apoptosis Transcription Factors and Regulators, Cancer Drug Targets, Tumor Suppressors

---

### Research Background

p53, also known as Tp53, is a DNA-binding protein which belongs to the p53 family. It contains transcription activation, DNA-binding, and oligomerization domains. p53 protein is expressed at low level in normal cells and at a high level in a variety of transformed cell lines, where it's believed to contribute to transformation and malignancy. p53 (TP53) is a transcription factor whose protein levels and post-translational modification state alter in response to cellular stress (such as DNA damage, hypoxia, spindle damage). Activation of p53 begins through a number of mechanisms including phosphorylation by ATM, ATR, Chk1 and MAPKs. MDM2 is a ubiquitin ligase that binds p53 and targets p53 for proteasomal degradation. Phosphorylation, p14ARF and USP7 prevent MDM2-p53 interactions, leading to an increase in stable p53 tetramers in the cytoplasm. Further modifications such as methylation and acetylation lead to an increase in Tp53 binding to gene specific response elements. Tp53 regulates a large number of genes (>100 genes) that control a number of key tumor suppressing functions such as cell cycle arrest, DNA repair, senescence and apoptosis. Whilst the activation of p53 often leads to apoptosis, p53 inactivation facilitates tumor progression. It is postulated to bind to a p53-binding site and activate expression of downstream genes that inhibit growth and/or invasion, and thus function as a tumor suppressor. Mutants of p53 that frequently occur in a number of different human cancers fail to bind the consensus DNA binding site, and hence cause the loss of tumor suppressor activity. Defects in TP53 are a cause of esophageal cancer, Li-Fraumeni syndrome, lung cancer and adrenocortical carcinoma. Cancer Immunotherapy/Immune Checkpoint/Immunotherapy/Targeted Therapy

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481