

Anti-CD32B/Fcgr2b Antibody (31868)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Cynomolgus
Conjugation:	Unconjugated
Clone:	31868
Purification:	Protein A

Applications

Verified Activity:	Immunochemical staining of cynomolgus CD32b in cynomolgus liver with mouse monoclonal antibody (1:60, formalin-fixed paraffin embedded sections).
Application:	IHC-P
Recommended	IHC-P: 1:50-1:200

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Cynomolgus CD32b / FCGR2B / Fc gamma RIIB Protein (TMPY-01971)
Antigen Species:	Cynomolgus
Synonyms:	Fc fragment of IgG, low affinity IIb, receptor (CD32)
Biology Area:	Fc Receptors, ITIM/ITAM Immunoreceptors and Related Molecules

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481