

Anti-CD32B/Fcgr2b Antibody (8R778)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Cynomolgus
Conjugation:	Unconjugated
Clone:	8R778
Purification:	Protein A

Applications

Verified Activity:	Flow cytometric analysis of Cynomolgus CD32 expression on Cynomolgus monocytes. Cells were stained with purified anti-CD32, then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable monocytes.
Application:	FCM
Recommended	FCM: 1:25-1:100

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Cynomolgus CD32b/FCGR2B/Fc gamma RIIB Protein (TMPY-01971)
Antigen Species:	Cynomolgus
Synonyms:	Fc fragment of IgG, low affinity IIb, receptor (CD32)
Biology Area:	Fc Receptors, ITIM/ITAM Immunoreceptors and Related Molecules

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481