

Anti-CD32B/Fcgr2b Antibody-FITC (7L654)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Cynomolgus
Conjugation:	FITC
Clone:	7L654
Purification:	Protein A

Applications

Verified Activity:	Flow cytometric analysis of Cynomolgus CD32 expression on Cynomolgus monocytes. Cells were stained with FITC-conjugated anti-CD32. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable monocytes.
Application:	FCM
Recommended	5 µl/Test, 0.1 mg/ml

Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Cynomolgus CD32b/FCGR2B/Fc gamma RIIB Protein (TMPY-01971)
Antigen Species:	Cynomolgus
Synonyms:	Fc fragment of IgG, low affinity IIb, receptor (CD32)
Biology Area:	Fc Receptors, ITIM/ITAM Immunoreceptors and Related Molecules

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481