

Anti-CD32A Antibody (6V174)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Cynomolgus
Conjugation:	Unconjugated
Clone:	6V174
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:1000-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Cynomolgus CD32a / Fc gamma RIIA protein (TMPY-01975)
Antigen Species:	Cynomolgus
Synonyms:	Fc fragment of IgG receptor IIa
Biology Area:	Fc Receptors, ITIM/ITAM Immunoreceptors and Related Molecules

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481