

## Anti-CD32A Antibody-FITC (3K614)

### Product Details

Ig Type:	Mouse IgG1
Reactivity:	Rhesus
Conjugation:	FITC
Clone:	3K614
Purification:	Protein A

### Applications

Verified Activity:	Flow cytometric analysis of rhesus FCGR2A(CD32A) expression on rhesus blood lymphocytes. Cells were stained with FITC-conjugated anti-rhesus FCGR2A(CD32A). The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of lymphocytes.
Application:	FCM
Recommended	5 µl/Test, 0.1 mg/ml

### Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Rhesus CD32a/FCGR2A/Fc gamma RIIA Protein (TMPY-02080)
Antigen Species:	Rhesus
Synonyms:	Fc fragment of IgG receptor IIa
Biology Area:	Fc Receptors, ITIM/ITAM Immunoreceptors and Related Molecules

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481