

## Anti-CD32A Antibody (5E468)

### Product Details

Ig Type:	Rabbit IgG
Reactivity:	Rhesus
Conjugation:	Unconjugated
Clone:	5E468
Purification:	Protein A

### Applications

Verified Activity:	Flow cytometric analysis of rhesus FCGR2A(CD32A) expression on rhesus spleen lymphocytes. Cells were stained with purified anti-rhesus FCGR2A(CD32A), then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.
Application:	FCM
Recommended	FCM: 1:25-1:100

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Rhesus CD32a/FCGR2A/Fc gamma RIIA Protein (TMPY-02080)
Antigen Species:	Rhesus
Synonyms:	Fc fragment of IgG receptor IIa
Biology Area:	Fc Receptors, ITIM/ITAM Immunoreceptors and Related Molecules

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481