

## Anti-CD32A Antibody (3M13)

### Product Details

Ig Type:	Rabbit IgG
Reactivity:	Rhesus
Conjugation:	Unconjugated
Clone:	3M13
Purification:	Protein A

### Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Rhesus CD32a / Fc gamma RIIA / FCGR2A Protein (TMPY-02080)
Antigen Species:	Rhesus
Synonyms:	Fc fragment of IgG receptor IIa
Biology Area:	Fc Receptors, ITIM/ITAM Immunoreceptors and Related Molecules

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481