

Anti-FCGRT & B2M Heterodimer Antibody (5Y16)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	5Y16
Purification:	Protein A

Applications

Verified Activity:	1. Immunochemical staining of human FcRn (FCGRT & B2M Heterodimer) in human skin with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections).
	2. Immunochemical staining of human FcRn (FCGRT & B2M Heterodimer) in human prostate with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections).
	3. Immunochemical staining of human FcRn (FCGRT & B2M Heterodimer) in human ovarian cancer with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections).
Application:	ELISA,IHC-P
Recommended	ELISA: 1:5000-1:10000; IHC-P: 1:100-1:500

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human FCGRT & B2M Heterodimer protein (TMPY-02082)
Antigen Species:	Human
Synonyms:	Fc fragment of IgG, receptor, transporter, alpha;Fc fragment of IgG, receptor, transporter, α
Biology Area:	Fc Receptors

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481