

Anti-IL-12 Antibody (3Q952)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3Q952
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human Human IL12A & IL12B Heterodimer protein (TMPY-02109)
Antigen Species:	Human
Synonyms:	IL12;IL12 p70;NFSK;NKSF;CLMF;interleukin 12A;P35;IMD29;IL-12B;NKSF1;NKSF2;IL-12A;CLMF2;IMD28

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481