

## Anti-IL-12 Antibody-HRP (9B977)

### Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	HRP
Clone:	9B977
Purification:	Protein A

### Applications

Application:	ELISA
Recommended	ELISA: 0.1-1 µg/ml

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human IL12 / IL12A & IL12B Protein (TMPY-02109)
Antigen Species:	Human
Synonyms:	interleukin 12A;IL-12B;NKSF1;NKSF2;IMD29;IL12 p70;NKSF;IL12;IL-12A;CLMF;NFSK;P35;IMD28;CLMF2

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481