

Anti-Integrin alpha 6 beta 1 Antibody (5J67)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	5J67
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human ITGA6 & ITGB1 Heterodimer protein (TMPY-02158)
Antigen Species:	Human
Synonyms:	integrin, $\alpha 6$;ITGA6B;Integrin $\alpha 6 \beta 1$;integrin, alpha 6;Integrin alpha 6 $\beta 1$;VLA-6;CD49f

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481