

Anti-Integrin alpha 5 beta 1 Antibody (5A186)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	5A186
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:1000-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human ITGA5 & ITGB1 Heterodimer protein (TMPY-02159)
Antigen Species:	Human
Synonyms:	FNRA;Integrin $\alpha 5 \beta 1$;VLA5A;VLA-5;integrin, $\alpha 5$ (fibronectin receptor, α polypeptide);Integrin alpha 5 $\beta 1$;CD49e;integrin, alpha 5 (fibronectin receptor, alpha polypeptide)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481