

Anti-CD3D & CD3E Heterodimer Antibody (2F74)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human, Cynomolgus
Conjugation:	Unconjugated
Clone:	2F74
Purification:	Protein A

Applications

Verified Activity:	<ol style="list-style-type: none">1. Immunochemical staining CD3E in human tonsil with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections).2. Immunochemical staining of human CD3E in cynomolgus macaque tonsil with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections).3. Immunochemical staining of human CD3E in cynomolgus macaque lymphnode with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections).4. Flow cytometric analysis of Human CD3E expression on human whole blood lymphocytes. Cells were stained with purified anti-Human CD3E, then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes.5. Immunochemical staining CD3 in cynomolgus spleen with rabbit monoclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.6. Immunochemical staining CD3 in cynomolgus lymph node with rabbit monoclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.
Application:	ELISA,ELISA(Det),FCM,IHC-P
Recommended	ELISA: 1:5000-1:10000; IHC-P: 1:100-1:500; FCM: 1:25-1:100; ELISA(Det): 1:1000-1:10000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human CD3D & CD3E Heterodimer protein (TMPY-02501)
Antigen Species:	Human
Synonyms:	IMD19;CD3d molecule, delta (CD3-TCR complex);T3D;CD3-DELTA;CD3d molecule, δ (CD3-TCR complex);CD3-δ

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481