

## Anti-CD3D & CD3E Heterodimer Antibody-APC (2F74)

### Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	APC
Clone:	2F74
Purification:	Protein A

### Applications

Verified Activity:	Flow cytometric analysis of Human CD3E expression on human whole blood lymphocytes. Cells were stained with APC-conjugated anti-Human CD3E. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes.
Application:	FCM
Recommended	5 µl/Test, 0.1 mg/ml

### Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human CD3D & CD3E Heterodimer protein (TMPY-02501)
Antigen Species:	Human
Synonyms:	T3D;CD3-δ;IMD19;CD3-DELTA;CD3d molecule, δ (CD3-TCR complex);CD3d molecule, delta (CD3-TCR complex)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481