

## Anti-CD3D & CD3E Heterodimer Antibody-PerCP (2F74)

### Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	PerCP
Clone:	2F74
Purification:	Protein A

### Applications

Verified Activity:	Flow cytometric analysis of Human CD3 expression on human peripheral blood lymphocytes. Cells were stained with PerCP-conjugated anti-Human CD3 and APC conjugated anti-Human CD19 (BD Pharmingen™ Cat. No. 555415). The dot plots were derived from events with the forward and side light-scatter characteristics of viable lymphocytes.
Application:	FCM
Recommended	10 µl/Test, 0.1 mg/ml

### Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human CD3D & CD3E Protein (TMPY-02501)
Antigen Species:	Human
Synonyms:	CD3d molecule, delta (CD3-TCR complex);IMD19;CD3-δ;T3D;CD3d molecule, δ (CD3-TCR complex);CD3-DELTA

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481