

Anti-CD3D & CD3E Heterodimer Antibody-HRP (2F74)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	HRP
Clone:	2F74
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 0.1-1 µg/ml

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human CD3D & CD3E Heterodimer protein (TMPY-02501)
Antigen Species:	Human
Synonyms:	CD3-δ;IMD19;CD3-DELTA;T3D;CD3d molecule, delta (CD3-TCR complex);CD3d molecule, δ (CD3-TCR complex)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481