

Anti-CD3D & CD3E Heterodimer Antibody-PE (2F74)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	PE
Clone:	2F74
Purification:	Protein A

Applications

Verified Activity:	Flow cytometric analysis of Human CD3E expression on human whole blood lymphocytes. Cells were stained with PE-conjugated anti-Human CD3E. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes.
Application:	FCM
Recommended	5 µl/Test, 0.1 mg/ml

Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human CD3D & CD3E Protein (TMPY-02501)
Antigen Species:	Human
Synonyms:	T3D;CD3-δ;CD3d molecule, δ (CD3-TCR complex);IMD19;CD3d molecule, delta (CD3-TCR complex);CD3-DELTA

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481