

Anti-CD3D & CD3E Heterodimer Antibody (7U199)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Mouse
Conjugation:	Unconjugated
Clone:	7U199
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Mouse CD3D & CD3E Heterodimer Protein (TMPY-03546)
Antigen Species:	Mouse
Synonyms:	IMD19;T3D;CD3d molecule, delta (CD3-TCR complex);CD3- δ ;CD3d molecule, δ (CD3-TCR complex);CD3-DELTA

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481