

## Anti-IL-23 Antibody (6L227)

### Product Details

Ig Type:	Mouse IgG2a
Reactivity:	Marmoset
Conjugation:	Unconjugated
Clone:	6L227
Purification:	Protein A

### Applications

Application:	ELISA
Recommended	ELISA: 1:1000-1:2000

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Marmoset IL-23 (IL23A & IL12B Heterodimer) Protein (TMPY-03689)
Antigen Species:	Marmoset
Synonyms:	IL23P19;IL-23;IL-23A;P19;SGRF;IL23

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481