


## Roridin L2

## Chemical Properties

CAS No. :	85124-22-7
Formula:	C <sub>29</sub> H <sub>38</sub> O <sub>9</sub>
Molecular Weight:	530.614
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Roridin L2 is a natural product for research related to life sciences. The catalog number is T41264 and the CAS number is 85124-22-7.
Targets(IC50)	Others,Endogenous Metabolite

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8846 mL	9.4231 mL	18.8462 mL
5 mM	0.3769 mL	1.8846 mL	3.7692 mL
10 mM	0.1885 mL	0.9423 mL	1.8846 mL
50 mM	0.0377 mL	0.1885 mL	0.3769 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481