

Anti-ASFV (strain Ba71V) p30 Antibody (2S577)

Product Details

Ig Type:	IgG1
Reactivity:	African swine fever virus
Conjugation:	Unconjugated
Clone:	2S577
Purification:	Protein G purified

Applications

Verified Activity:	The Binding Activity of African swine fever virus Phosphoprotein p30(Ba71V-93). Activity: Measured by its binding ability in a functional ELISA. Immobilized African swine fever virus Phosphoprotein p30(Ba71V-93) protein at 2 µg/ml can bind African swine fever virus Phosphoprotein p30(Ba71V-93) Antibody, the EC50 is 0.01982 to 0.02784 ug/ml.
Application:	ELISA

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: PhosphoProtein p30 Protein (1-204AA)
Antigen Species:	African swine fever virus
Gene ID:	22220322
Uniprot ID:	P34204
Biology Area:	Others

Research Background

Modifies the subcellular distribution of heterogeneous nuclear ribonucleoprotein K (HNRNPK) and may contribute to modulate HNRNPK functions related to processing and export of mRNAs during ASFV infection.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481