

Anti-Fetuin A Antibody (7B541)

Product Details

| | |
|---------------|--------------------|
| Ig Type: | IgG1 |
| Reactivity: | Human |
| Conjugation: | Unconjugated |
| Clone: | 7B541 |
| Purification: | Protein G purified |

Applications

| | |
|--------------------|---|
| Verified Activity: | 1. Immunohistochemical of paraffin-embedded human liver cancer using TMAH-00043 at dilution of 1:200 |
| | 2. Immunohistochemical of paraffin-embedded human breast cancer using TMAH-00043 at dilution of 1:200 |
| Application: | ELISA,IHC |

Properties

| | |
|----------------------|---|
| Purity: | >95% |
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------------|--|
| Immunogen: | Recombinant Protein: Human AHSG |
| Antigen Species: | Human |
| Gene ID: | 197 |
| Uniprot ID: | P02765 |
| Synonyms: | FETUA;HSGA; α -2-HS-glycoprotein;AHS;A2HS;alpha-2-HS-glycoprotein |
| Biology Area: | Cardiovascular |

Research Background

Promotes endocytosis, possesses opsonic properties and influences the mineral phase of bone. Shows affinity for calcium and barium ions.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481