

Anti-AMBP Antibody (8J942)

Product Details

Ig Type:	IgG3
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	8J942
Purification:	Protein G purified

Applications

Verified Activity:	Immunohistochemical of paraffin-embedded human liver tissue using TMAH-00058 at dilution of 1:200
Application:	ELISA,IHC

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human Protein AMBP
Antigen Species:	Human
Uniprot ID:	P02760
Synonyms:	IATIL;HI30;HCP;ITI;A1M;EDC1;alpha-1-microglobulin/bikunin precursor;UTI;α-1-microglobulin/bikunin precursor;ITILC;ITIL
Biology Area:	Cardiovascular

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481