

Anti-B2M/beta 2-Microglobulin Antibody (8H120)

Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	8H120
Purification:	Protein G purified

Applications

Verified Activity:	Immunohistochemical of paraffin-embedded human tonsil tissue using TMAH-00100 at dilution of 1:200
Application:	ELISA, IHC

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: human Beta-2-microglobulin Protein
Antigen Species:	Human
Uniprot ID:	P61769
Synonyms:	β -2-microglobulin;Beta-2-microglobulin;B2M;B2M/ β 2-Microglobulin
Biology Area:	Cardiovascular

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481