

Anti-CALCA Antibody (4W385)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	4W385
Purification:	Protein G purified

Applications

Verified Activity:	Immunohistochemical of paraffin-embedded human thyroid tissues using TMAH-00133 at dilution of 1:200
Application:	ELISA, IHC

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human PCT
Antigen Species:	Human
Gene ID:	796
Uniprot ID:	P01258
Synonyms:	calcitonin-related polypeptide α ;CGRP;CGRP-I;CGRP1;CT;Procalcitonin;KC;calcitonin-related polypeptide alpha;PCT;CALC1
Biology Area:	Neuroscience

Research Background

Calcitonin causes a rapid but short-lived drop in the level of calcium and phosphate in blood by promoting the incorporation of those ions in the bones.; Katalcalcin is a potent plasma calcium-lowering peptide.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481