

## Anti-CCL2 Antibody (5R247)

### Product Details

Ig Type:	IgG2b
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	5R247
Purification:	Protein G purified

### Applications

Verified Activity:	Immunohistochemical of paraffin-embedded human pancreas using TMAH-00150 at dilution of 1:50
Application:	ELISA, IHC

### Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human C-C motif chemokine 2
Antigen Species:	Human
Uniprot ID:	P13500
Synonyms:	SCYA2;C-C motif chemokine 2;HSMCR30;MCP-1;MCAF;SMC-CF;HC11;GDGF-2;CCL2;MCP1
Biology Area:	Immunology

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481