

Anti-CCR5 Antibody (6D513)

Product Details

Ig Type:	hIgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6D513
Purification:	Affinity-chromatography

Applications

Verified Activity:	<ol style="list-style-type: none">1. The Binding Activity of Human CCR5 with Anti-CCR5 Recombinant Antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CCR5 at 10 µg/mL can bind Anti-CCR5 recombinant antibody. The EC50 is 1.099-1.287 ng/mL. The VLPs is negative control.2. The purity of CCR5 was greater than 95% as determined by SEC-HPLC
Application:	ELISA, FCM
Recommended	FCM:1:50-1:200.

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human CCR5 Protein
Antigen Species:	Human
Uniprot ID:	P51681
Synonyms:	CCCKR5;CKR-5;IDDM22;CMKBR5;CD195;CKR5;CC-CKR-5;CCR-5
Biology Area:	Neuroscience, Cancer, Cardiovascular, Immunology, Microbiology, Signal transduction

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481