

## Anti-CCR8 Antibody (3G66)

## Product Details

Ig Type:	Human IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3G66
Purification:	Affinity-chromatography

## Applications

Verified Activity:	Untransfected HEK293T cells surface (green line) and transfected Human CCR8 HEK293T stable cells surface (red line) were stained with anti-CCR8 antibody (2 $\mu$ g/1*10 <sup>6</sup> cells), washed and then followed by FITC-conjugated anti-Human IgG Fc antibody and analyzed with flow cytometry.
Application:	FCM

## Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	Recombinant Protein: Human CCR8 Protein
Antigen Species:	Human
Gene ID:	1237
Uniprot ID:	P51685
Synonyms:	CCR-8;TER-1;CMKBRL2;GPR-CY6;CKRL1;CD198;CKR-L1;GPRCY6;TER1;CCR8;CY6;CDw198; ChemR1;CMKBR8
Biology Area:	Immunology, Microbiology

## Research Background

Receptor for the chemokine CCL1/SCYA1/I-309. May regulate monocyte chemotaxis and thymic cell line apoptosis. Alternative coreceptor with CD4 for HIV-1 infection.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481