

Anti-CCR8 Antibody (6N151)

Product Details

Ig Type:	Human IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6N151
Purification:	Affinity-chromatography

Applications

Verified Activity:	1. Untransfected CHO-K1*10 ⁶ cells, washed and then followed by FITC-conjugated anti-Human IgG Fc antibody and analyzed with flow cytometry. 2. The Binding Activity of CCR8 with Anti-CCR8 recombinant antibody. Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CCR8 at 5 µg/mL can bind Anti-CCR8 recombinant antibody, the EC50 is 11.20-15.63 ng/mL.
Application:	FCM

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human CCR8 Protein
Antigen Species:	Human
Gene ID:	1237
Uniprot ID:	P51685
Synonyms:	CY6;ChemR1;CDw198;TER1;CMKBRL2;CD198;TER-1;CKR-L1;CCR8;CKRL1;GPR-CY6;CCR-8;CMKBR8;GPRCY6
Biology Area:	Immunology

Research Background

Receptor for the chemokine CCL1/SCYA1/I-309. May regulate monocyte chemotaxis and thymic cell line apoptosis. Alternative coreceptor with CD4 for HIV-1 infection.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481