

## Anti-CCR8 Antibody (6H752)

### Product Details

Ig Type:	Human IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6H752
Purification:	Affinity-chromatography

### Applications

Verified Activity:	Untransfected CHO-K1*10 <sup>6</sup> cells, washed and then followed by FITC-conjugated anti-Human IgG Fc antibody and analyzed with flow cytometry.
Application:	FCM

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human CCR8 Protein
Antigen Species:	Human
Gene ID:	1237
Uniprot ID:	P51685
Synonyms:	CKR-L1;GPRCY6;TER1;CY6;CMKBRL2;CCR8;CKRL1;CD198;CDw198;CCR-8;ChemR1;TER-1;CMKBR8;GPR-CY6
Biology Area:	Immunology

### Research Background

Receptor for the chemokine CCL1/SCYA1/I-309. May regulate monocyte chemotaxis and thymic cell line apoptosis. Alternative coreceptor with CD4 for HIV-1 infection.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481