

Anti-CD19 Antibody (2A873)

Product Details

Ig Type:	Mouse IgG2a
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	2A873
Purification:	Affinity-chromatography

Applications

	Western Blot
	-Positive WB detected in: Raji whole cell lysate, JK whole cell lysate
Verified Activity:	-All lanes: CD19 antibody at 1:1000
	-Secondary: Goat polyclonal to mouse IgG at 1/50000 dilution
	-Predicted band size: 62 kDa
	-Observed band size: 62 kDa
Application:	ELISA, WB
Recommended	WB: 1:500-1:5000.

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human CD19 Protein
Antigen Species:	Human
Gene ID:	930
Uniprot ID:	P15391
Synonyms:	CVID3; B4; CD19 molecule
Biology Area:	Immunology, Stem cells

Research Background

Functions as coreceptor for the B-cell antigen receptor complex (BCR) on B-lymphocytes. Decreases the threshold for activation of downstream signaling pathways and for triggering B-cell responses to antigens. Activates signaling pathways that lead to the activation of phosphatidylinositol 3-kinase and the mobilization of intracellular Ca(2+) stores. Is not required for early steps during B cell differentiation in the blood marrow. Required for normal differentiation of B-1 cells. Required for normal B cell differentiation and proliferation in response to antigen challenges. Required for normal levels of serum immunoglobulins, and for production of high-affinity antibodies in response to antigen challenge.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481