

Anti-CDC45 Antibody (8Z278)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	8Z278
Purification:	Affinity-chromatography

Applications

Application:	ELISA
--------------	-------

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthetic peptide: Human CDC45
Antigen Species:	Human
Gene ID:	8318
Uniprot ID:	O75419
Synonyms:	Cell division control protein 45 homolog; CDC45 related protein; CDC45 like; CDC45 cell division cycle 45 like; cell division cycle 45 (S.cerevisiae homolog) like; CDC 45L; CDC45L2; CDC45; PORC PI 1; PORC-PI-1; Cell division cycle 45 like; CDC 45; Cell division cycle 45 like 2
Biology Area:	Epigenetics and Nuclear Signaling

Research Background

Required for initiation of chromosomal DNA replication.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481