

Anti-CDH1 Antibody (7T289)

Product Details

Ig Type:	hIgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7T289
Purification:	Affinity-chromatography

Applications

Verified Activity:	The Binding Activity of Human CDH1 with Anti-CDH1 recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CDH1 at 2 µg/mL can bind Anti-CDH1 recombinant antibody.The EC50 is 5.151-6.358 ng/mL.
Application:	ELISA
Recommended	ELISA:1:5000-1:50000.

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human CDH1 Protein
Antigen Species:	Human
Uniprot ID:	P12830
Synonyms:	Arc-1;E-cad;UVO;CDH1;E-Cadherin;ECAD;CDHE;cadherin 1, type 1, E-cadherin (epithelial); CD324;LCAM
Biology Area:	Cancer, Developmental biology, Signal transduction

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481