

Anti-CEACAM5 Antibody (6V639)

Product Details

Ig Type:	hIgG1
Reactivity:	Cynomolgus
Conjugation:	Unconjugated
Clone:	6V639
Purification:	Affinity-chromatography

Applications

Verified Activity:	The Binding Activity of Cynomolgus CEACAM5 with Anti-CEACAM5 Recombinant Antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Cynomolgus CEACAM5 at 2µg/mL can bind Anti-CEACAM5 recombinant antibody. The EC50 is 3.572-4.044 ng/mL.
Application:	ELISA
Recommended	ELISA:1:5000-1:50000.

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human CEACAM5 Protein
Antigen Species:	Human
Uniprot ID:	P06731
Synonyms:	carcinoembryonic antigen-related cell adhesion molecule 5;CD66e;CEA
Biology Area:	Cancer, Tags & Cell Markers, Immunology, Stem cells

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481