

## Anti-CEACAM8 Antibody (9F528)

### Product Details

Ig Type:	Human IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	9F528
Purification:	Affinity-chromatography

### Applications

Verified Activity:	The Binding Activity of Human CEACAM8 with Anti-CEACAM8 Recombinant Antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CEACAM8 at 2 µg/mL can bind Anti-CEACAM8 recombinant antibody, the EC50 is 19.38-21.68 ng/mL.
Application:	ELISA

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human CEACAM8 Protein
Antigen Species:	Human
Gene ID:	1088
Uniprot ID:	P31997
Biology Area:	Immunology

### Research Background

Cell surface glycoprotein that plays a role in cell adhesion in a calcium-independent manner. Mediates heterophilic cell adhesion with other carcinoembryonic antigen-related cell adhesion molecules, such as CEACAM6. Heterophilic interaction with CEACAM8 occurs in activated neutrophils.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481