

## Anti-Chenodeoxycholic acid Antibody (2K943)

### Product Details

|               |                    |
|---------------|--------------------|
| Ig Type:      | IgG2b              |
| Conjugation:  | Unconjugated       |
| Clone:        | 2K943              |
| Purification: | Protein G purified |

### Applications

|              |       |
|--------------|-------|
| Application: | ELISA |
|--------------|-------|

### Properties

|                      |   |
|----------------------|---|
| Purity:              | >95%  |
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping:            | Shipping with blue ice.   |

### Antigen Details

|               |   |
|---------------|---|
| Immunogen:    | Chenodeoxycholic acid-BSA conjugated  |
| Synonyms:     | chenodesoxycholic acid;chenix;CDCA;chenodeoxycholate;chenic acid;chenofalk;<br>gallodesoxycholic acid;chendol;Chenodiol |
| Biology Area: | Others  |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481