

Anti-CHGA Antibody (6C108)

Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6C108
Purification:	Protein G purified

Applications

1. Immunohistochemical of paraffin-embedded Human small intestine tissue using TMAH-00270 at dilution of 1:200.
2. Immunohistochemical of paraffin-embedded Human pancreas tissue using TMAH-00270 at dilution of 1:200.
3. Western Blot

Verified Activity:

- All lanes: Chromogranin-A transfected E.coli lysate
- Lane 1: Mouse Anti--6*his monoclonal antibody at 1µg/ml
- Lane 2 : Mouse anti- Chromogranin-A monoclonal antibody at 1µg/ml
- Predicted band size : 52Kd
- Observed band size : 52kd

Application: ELISA,IHC,WB

Properties

Purity: >95%

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: CHGA Protein
Antigen Species:	Human
Uniprot ID:	P10645
Biology Area:	Neuroscience

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481