

## Anti-Clenbuterol Antibody (8E543)

### Product Details

Ig Type: IgG1  
Conjugation: Unconjugated  
Clone: 8E543  
Purification: Protein G purified

### Applications

Application: ELISA

### Properties

Purity: >95%  
Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

### Antigen Details

Immunogen: Clenbuterol-BSA conjugated  
Biology Area: Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481