

Anti-CST3 Antibody (1W125)

Product Details

Ig Type:	IgG2b
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	1W125
Purification:	Protein G purified

Applications

Verified Activity:	1. Immunohistochemical of paraffin-embedded human prostate tissue using TMAH-00307 at dilution of 1:200
	2. Immunohistochemical of paraffin-embedded human kidney tissue using TMAH-00307 at dilution of 1:200
Application:	ELISA, IHC

Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human Cystatin C Protein
Antigen Species:	Human
Uniprot ID:	P01034
Synonyms:	Cystatin-3;CST3;Post-gamma-globulin;Neuroendocrine basic polypeptide;Gamma-trace;Post- γ -globulin; γ -trace;Cystatin-C
Biology Area:	Tags & Cell Markers

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481