

## Anti-CTGF/CCN2 Antibody (7D485)

### Product Details

|               |                         |
|---------------|-------------------------|
| Ig Type:      | hIgG1                   |
| Reactivity:   | Macaca fascicularis     |
| Conjugation:  | Unconjugated            |
| Clone:        | 7D485                   |
| Purification: | Affinity-chromatography |

### Applications

|                    |   |
|--------------------|---|
| Verified Activity: | The Binding Activity of Rhesus macaque CTGF with Anti-CTGF Recombinant Antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Rhesus macaque CTGF at 1 µg/mL can bind Anti-CTGF recombinant antibody. The EC50 is 0.6886-0.7825 ng/mL. |
| Application:       | ELISA   |
| Recommended        | ELISA:1:5000-1:50000.   |

### Properties

|                      |   |
|----------------------|---|
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping:            | Shipping with blue ice.   |

### Antigen Details

|                  |   |
|------------------|---|
| Immunogen:       | Recombinant Protein: Human CTGF Protein |
| Antigen Species: | Human                                   |
| Uniprot ID:      | P29279                                  |
| Biology Area:    | Cardiovascular, Immunology              |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481