

Anti-CTLA4 Antibody (3K370)

Product Details

Ig Type:	hIgG1
Reactivity:	Peptostreptococcus magnus
Conjugation:	Unconjugated
Clone:	3K370
Purification:	Affinity-chromatography

Applications

Verified Activity:	1. The Binding Activity of Peptostreptococcus magnus protein L with anti-CTLA4 Recombinant Antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Peptostreptococcus magnus protein L at 2 µg/mL can bind anti-CTLA4 recombinant antibody. The EC50 is 1.764-2.413 ng/mL.
	2. The Binding Activity of Peptostreptococcus magnus protein L with Anti-CTLA4 Recombinant Antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Peptostreptococcus magnus protein L at 2 µg/mL can bind Anti-CTLA4 recombinant antibody. The EC50 is 1.601-1.944 ng/mL.
Application:	ELISA
Recommended	ELISA:1:5000-1:50000.

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human proetin L Protein
Antigen Species:	Human
Uniprot ID:	P16410
Synonyms:	CTLA4;ALP55;IDDM12;CELIAC3;CD152;ICOS;GRD4;GSE
Biology Area:	Immunology, Stem cells

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481