

## Anti-Dehydroepiandrosterone(DHEA) Antibody (4I201)

### Product Details

Ig Type:	IgG1
Conjugation:	Unconjugated
Clone:	4I201
Purification:	Protein G purified

### Applications

Application:	ELISA
--------------	-------

### Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Dehydroepiandrosterone-BSA conjugated
Synonyms:	17-hormoforin;diandrone;psicosterone;trans-dehydroandrosterone;DHEA;diandron;androstenolone;dehydroisoandrosterone;prasterone
Biology Area:	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481