

Anti-DM1 Antibody (2M450)

Product Details

Ig Type:	IgG2a
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	2M450
Purification:	Affinity-chromatography

Applications

Verified Activity:	The Binding Activity of Trastuzumab emtansine (Ado-Trastuzumab emtansine) and Cantuzumab ravtansine (Synonyms:IMGN242) with Anti-DM1 antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Trastuzumab emtansine (Ado-Trastuzumab emtansine) and Cantuzumab ravtansine (Synonyms:IMGN242) at 2 µg/mL can bind Anti-DM1 antibody, the EC50 is 2.697 to 2.984 ng/mL and no.
Application:	ELISA

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	DM1-KLH
Synonyms:	Mertansine;DM 1;Maytansinoid DM1
Biology Area:	DM1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481