

Anti-DPEP3 Antibody (2P258)

Product Details

Ig Type:	mIgG1
Reactivity:	Human, Macaca fascicularis
Conjugation:	Unconjugated
Clone:	2P258
Purification:	Affinity-chromatography

Applications

Verified Activity:	1. The Binding Activity of Human DPEP3 with Anti-DPEP3 recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human DPEP3 at 2 µg/mL can bind Anti-DPEP3 recombinant antibody. The EC50 is 6.841-7.498 ng/mL. 2. The Binding Activity of Macaca fascicularis DPEP3 with Anti-DPEP3 recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Macaca fascicularis DPEP3 at 2 µg/mL can bind Anti-DPEP3 recombinant antibody. The EC50 is 7.817-8.936 ng/mL.
Application:	ELISA
Recommended	ELISA:1:5000-1:50000.

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human DPEP3 Protein
Antigen Species:	Human
Uniprot ID:	Q9H4B8
Synonyms:	DPEP 3;UNQ834/PRO1772;Dipeptidase 3
Biology Area:	Cell Biology

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481