

## Anti-DXD Antibody-HRP (9W151)

## Product Details

Ig Type:	IgG2a
Conjugation:	HRP
Clone:	9W151
Purification:	Protein G purified

## Applications

Verified Activity:	<p>1. The Binding Activity of ADC-1(DXD) with Anti-DXD-HRP antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized ADC-1(DXD) at 2 µg/mL can bind Anti-DXD-HRP antibody, the EC50 is 8.653 to 12.39 ng/mL.</p> <p>2. The Binding Activity of ADC-2(DXD) with Anti-DXD-HRP antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized ADC-2(DXD) at 2 µg/mL can bind Anti-DXD-HRP antibody, the EC50 is 9.854 to 11.57 ng/mL.</p> <p>3. The Binding Activity of T-DXd(DS-8201) with Anti-DXD-HRP antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized T-DXd(DS-8201) at 2 µg/mL can bind Anti-DXD-HRP antibody, the EC50 is 6.726 to 8.530 ng/mL.</p>
Application:	ELISA

## Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	DXD-BSA conjugated
Synonyms:	Exatecan derivative
Biology Area:	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481