

## Anti-EDAR Antibody (6J576)

### Product Details

Ig Type:	Mouse IgG2a
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6J576
Purification:	Affinity-chromatography

### Applications

Application:	ELISA, IHC
Recommended	IHC:1:20-1:200.

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human EDAR Protein
Antigen Species:	Human
Gene ID:	#N/A
Uniprot ID:	Q9UNE0
Synonyms:	EDA1R;ED5;ECTD10A;ectodysplasin A receptor;ED3;EDA-A1R;ECTD10B;ED1R;EDA3;HRM1;DL
Biology Area:	Epigenetics and Nuclear Signaling, Signal transduction

### Research Background

#N/A

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481