

Anti-Eribulin Antibody (8T987)

Product Details

Ig Type:	IgG2a
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	8T987
Purification:	Affinity-chromatography

Applications

Verified Activity:	1. The Binding Activity of Farletuzumab ecteribulin (Synonyms: MORAb-202) with Anti-Eribulin antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Farletuzumab ecteribulin (Synonyms: MORAb-202) at 2 µg/mL can bind Anti-Eribulin antibody, the EC50 is 0.8548 to 1.015 ng/mL. 2. The purity of Eribulin was greater than 95% as determined by SEC-HPLC
Application:	ELISA
Recommended	ELISA:1:5000-1:50000.

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Eribulin-BSA
Synonyms:	Eribulin;B1939;ER-086526;E7389
Biology Area:	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481