

## Anti-FABP3 Antibody (7D420)

## Product Details

Ig Type:	IgG2a
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7D420
Purification:	Protein G purified

## Applications

1. Immunohistochemical of paraffin-embedded human breast cancer tissue using TMAH-00416 at dilution of 1:200
2. Immunohistochemical of paraffin-embedded human kidney tissue using TMAH-00416 at dilution of 1:200
3. Western Blot

Verified Activity:	-All lanes: Mouse anti- Human Fatty acid-binding protein monoclonal antibody at 1µg/ml -Lane 1:mouse heart muscle lysate -Lane 2:rat heart muscle lysate -Lane 3:FABP3 transfected SF9 cell lysate -Predicted band size : 15kd -Observed band size : 15kd
--------------------	--

Application: ELISA, WB, IHC

## Properties

Purity:	>95%
Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	Recombinant Protein: Human Fatty acid-binding Protein
Antigen Species:	Human
Uniprot ID:	P05413
Synonyms:	Mammary-Derived Growth Inhibitor;MDGIMuscle Fatty Acid-Binding Protein;Fatty Acid-Binding Protein 3;Fatty Acid-Binding Protein Heart;MDGI;Heart-Type Fatty Acid-Binding Protein;H-FABP; FABP3;FABP11;M-FABP
Biology Area:	Cardiovascular

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481